



970-686-2266 | WWW.QUORUMPROSTHETICS.COM

SCANNING

All 3D printed sockets start with a scan of either the residual limb, a plaster cast or a well fitting socket.



Click Dial
OR

Icarus Push Button



QUATRO PIN-LOCK



PATENT NO.

10004614B1, 10406003B2, 11844667B2

CONFIDENCE & STABILITY

The Quatro Pin-Lock is designed to work in tandem with our comfort cells. These cells can be easily interchanged on the socket rails. This allows the user to quickly change cells to maximize stability and comfort.



Comfort Cells
PATENT PENDING

DESIGN

Our highly skilled team of engineers will use the digital scan to make the socket 100% custom for you.

Quatro sockets are rigorously tested to ensure you will have a lightweight and incredibly strong socket for many years.

The Quatro Pin-Lock is designed with four zones of proportional compression control.

Carbon Reinforced

Above Knee & Below Knee sockets are available now!

MADE IN THE USA

Each 3D printed device is manufactured in the United States of America using the highest quality materials.